

InVision-Plus® with Neutral Advantage Technology

DIRECTIONS - Use aseptic technique Wash hands thoroughly with bactericidal soap before each procedure. Use gloves if required by your Healthcare provider.

Inspect device. Discard device if end caps are missing or loose in package. Remove protective coverings as assembly progresses.

InVision-Plus® with Neutral Advantage Technology

The InVision-Plus® is accessed by a standard male-luer syringe or I.V. set connector. It is recommended that a male luer-lock connector be used with the InVision-Plus® for a secure connection. The InVision-Plus® is bi-directional and luer-locking. The InVision-Plus® will automatically close when the syringe or male-luer I.V. set connector is removed. The InVision-Plus® has been tested with low pressure power injectors having a maximum pressure of 325 psi and a maximum flow rate of 10mL/sec. Connect only with other devices rated for high pressure with a luer lock connection when performing power injections.

InVision-Plus® with Neutral Advantage Technology

1. Peel Package: Peel back the top of the Peel Package. Remove end cap. Do not touch the tip of the InVision-Plus®. Attach the InVision-Plus® to the desired vascular access device by rotating the InVision-Plus® in a clockwise direction until secure.

"Touch-Free" Container Package: Peel back the lid on the "Touch-Free" Package Container. Do not touch the tip of the InVision-Plus®. Attach the InVision-Plus® to the desired vascular access device by rotating the "Touch-Free" Package Container in a clockwise direction until secure. Pull back the "Touch-Free" Container until clear of the port and dispose. Note: device does not have end caps

2. To access the InVision-Plus® swab the septum and thread area with a facility approved disinfectant.
3. Position the syringe or I.V. set male-luer connector tip on the InVision-Plus® Septum.
4. Push forward, rotate syringe or I.V. set male luer-lock connector clockwise until secure.
5. To disconnect the syringe or I.V. set male-luer connector, simply rotate counterclockwise. The InVision-Plus® will automatically close.
6. The InVision-Plus®/ InVision-Plus® RED® has 0.027mL priming volume.
The InVision-Plus® Junior®/ InVision-Plus® EPI™ has 0.022mL priming volume.

InVision-Plus® CATHETER EXTENSION SETS

Prime Set

1. Open clamp, prime to expel air from set.
2. Close clamp
3. Attach set to vascular access device.
4. Ensure that all connections are tightened and secure before use.









NOTES

Flush InVision-Plus® after each use per institutional protocol using normal saline, or normal saline and heparin. It is recommended that this device be changed per CDC guidelines or per validated facility protocol. The CDC recommends that tubing used to administer lipid emulsions and TPN solutions should be changed every 24 hours. The IV connector and extension set should also be changed at this time. Observe appropriate infection control procedures. The InVision-Plus® EPI™ is yellow in color by design as a visual identifier to be used on epidural catheters.

Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTIONS

The InVision-Plus® may not be compatible with every male-luer connector. Some I.V. sets/syringes utilize a male-luer connector with a lumen opening diameter too small to be used on the InVision-Plus® which may cause damage to the InVision-Plus® resulting in septum damage, reduction or loss of fluid flow. Do not use any male-luer device on the InVision-Plus®/ InVision-Plus® RED® unless certain that the male-luer connector opening is 0.064" or larger. Do not use any male-luer device on the InVision-Plus® Junior®/ InVision-Plus® EPI™ unless certain that the male-luer connector opening is 0.057" or larger. Do not use needles with connector. Do not use luer lock end caps on connector.

 Radiation Sterilized  Ethylene Oxide Sterilized  Do Not Reuse  Do Not Re-Sterilize
 DEHP Free  Latex Free  See Directions for Use  Non Pyrogenic Fluidpathway RX Only

Product Information: (512) 301-7334 www.rymedtech.com

InVision-Plus®, RED®, Junior® and EPI™ are trademarks of RyMed Technologies, LLC Patents: #6,113,068 - #6,299,131 - #6,994,315 - #7,530,546B2 - #8,096,525B2 and other patents pending © RyMed Technologies, LLC 2014 3RT506N