

# InVision-Plus®

## Product Catalog

---



**InVision-Plus®**  
*with Neutral Advantage Technology*  
*Intraluminal Protection Systems*

## InVision-Plus® with Neutral Advantage Technology Intraluminal Protection Systems

Designed for I.V. therapy procedures, blood administration and blood sampling.  
All products and packaging are not made with natural rubber latex.

### Catalog/Part Number      Product Description


**RYM – 2001**      **InVision-Plus® RED®**  
PEEL POUCH - Used for sterile field procedures  
Intended for Visual Identification  
PV = 0.027mL



**RYM – 5000**      **InVision-Plus®**  
TOUCH-FREE™ Container Package  
Intended for Use on All Patients  
PV = 0.027mL



**RYM – 5000CL**      **InVision-Plus® Clear**  
TOUCH-FREE™ Container Package  
Intended for Use on All Patients  
PV = 0.027mL




**RYM – 5001**      **InVision-Plus®**  
PEEL POUCH - Used for sterile field procedures  
Intended for Use on All Patients  
PV = 0.027mL



**RYM – 5001CL**      **InVision-Plus® Clear**  
PEEL POUCH - Used for sterile field procedures  
Intended for Use on All Patients  
PV = 0.027mL



**RYM – 7000**      **InVision-Plus CS®**  
**w/Chlorhexedine+Silver Engineering**  
TOUCH-FREE™ Container Package  
Intended for Use on All Patients  
PV = 0.027mL



**RYM – 7001**      **InVision-Plus CS®**  
**w/Chlorhexedine+Silver Engineering**  
PEEL POUCH - Used for sterile field procedures  
Intended for Use on All Patients  
PV = 0.027mL



**RYM – 8001CL**      **InVision-Plus® Junior® Clear**  
PEEL POUCH - Used for sterile field procedures  
**Intended for Use on Neonatal & Pediatric Patients**  
PV = 0.022mL









Packaging: 200 Units per Box - 5 Boxes per Case - TOTAL: 1,000 UNITS PER CASE

## Catheter Extension Sets

### With InVision-Plus® with Neutral Advantage Technology Intraluminal Protection System

NOTE: All catheter extension sets are made of soft, kink-resistant, NON-DEHP PVC Tubing.  
All products and packaging are not made with natural rubber latex.

Catalog/Part Number	Product Description	
<b>RYM – 5104B</b>	<b>4" Microbore Bifurcated Catheter Extension Set</b> <b>(2) Pre-attached InVision-Plus®</b> (2) Slide clamps, rotating male-luer lock Note: InVision-Plus® is not bonded to set APV = 0.41mL	
<b>RYM – 5104T</b>	<b>4" Microbore T-Port Catheter Extension Set</b> <b>(1) Pre-attached InVision-Plus® Clear</b> Slide clamp, rotating male-luer lock APV = 0.26mL	
<b>RYM – 5107</b>	<b>7" Microbore Catheter Extension Set</b> <b>(1) Pre-attached InVision-Plus®</b> Slide clamp, fixed male-luer lock Note: InVision-Plus® is not bonded to set APV = 0.27mL	
<b>RYM – 5107CL</b>	<b>7" Microbore Catheter Extension Set</b> <b>(1) Pre-attached InVision-Plus® Clear</b> Slide clamp, fixed male-luer lock Note: InVision-Plus® is not bonded to set APV = 0.27mL	
<b>RYM – 5107HPB</b>	<b>7" High Pressure Microbore Bifurcated Cath. Ext. Set*</b> <b>(2) Pre-attached InVision-Plus®</b> *High Pressure Catheter Extension Set rated 10mL/sec w/max pressure of 325 psi (2) Slide clamps, rotating male-luer lock Note: InVision-Plus® is not bonded to set APV = 0.45mL	
<b>RYM – 5107HPBCL</b>	<b>7" High Pressure Microbore Bifurcated Cath. Ext. Set*</b> <b>(2) Pre-attached InVision-Plus® Clear</b> *High Pressure Catheter Extension Set rated 10mL/sec w/max pressure of 325 psi (2) Slide clamps, rotating male-luer lock Note: InVision-Plus® is not bonded to set APV = 0.45mL	

Packaging: 100 Units per Box - 4 Boxes per Case - TOTAL: 400 UNITS PER CASE

## Catheter Extension Sets


With InVision-Plus® with Neutral Advantage Technology

*Intraluminal Protection System*


NOTE: All catheter extension sets are made of soft, kink-resistant, NON-DEHP PVC Tubing.  
All products and packaging are not made with natural rubber latex.

Catalog/Part Number	Product Description
---------------------	---------------------


<b>RYM – 5307HPCL</b>	<p><b>7” High Pressure Macrobore Catheter Extension Set*</b>  <b>(1) Pre-attached InVision-Plus® Clear</b>                      *High Pressure Catheter Extension Set rated 10mL/sec w/max pressure of 325 psi                      Slide clamp, universal rotating male-luer lock                      Note: InVision-Plus® is not bonded to set                      APV = 0.42mL</p>
-----------------------	---




<b>RYM – 5307HPU</b>	<p><b>7” High Pressure Macrobore Catheter Extension Set*</b>  <b>(1) Pre-attached InVision-Plus®</b>                      *High Pressure Catheter Extension Set rated 10mL/sec w/max pressure of 325 psi                      Slide clamp, universal rotating male-luer lock                      Note: InVision-Plus® is not bonded to set                      APV = 0.42mL</p>
----------------------	---




<b>RYM – 5315</b>	<p><b>15” Macrobore Catheter Extension Set</b>  <b>(1) Pre-attached InVision-Plus®</b>                      Slide clamp, fixed male-luer lock                      Note: InVision-Plus® is not bonded to set                      APV = 0.63mL</p>
-------------------	--




<b>RYM – 5407</b>	<p><b>7” Standardbore Catheter Extension Set</b>  <b>(1) Pre-attached InVision-Plus®</b>                      Slide clamp, fixed male-luer lock                      Note: InVision-Plus® is not bonded to set                      APV = 1.32mL</p>
-------------------	--



<b>RYM – 5407CL</b>	<p><b>7” Standardbore Catheter Extension Set</b>  <b>(1) Pre-attached InVision-Plus® Clear</b>                      Slide clamp, fixed male-luer lock                      Note: InVision-Plus® is not bonded to set                      APV = 1.32mL</p>
---------------------	--



<b>RYM – 5412Y</b>	<p><b>12” Standardbore IV Extension Set</b>  <b>(1) InVision-Plus® Clear Y-Injection Site</b>                      Slide clamp, fixed male-luer lock                      APV = 2.26mL</p>
--------------------	--







Packaging: 100 Units per Box - 4 Boxes per Case - TOTAL: 400 UNITS PER CASE

## Catheter Extension Sets







### With InVision-Plus® with Neutral Advantage Technology Intraluminal Protection System

NOTE: All catheter extension sets are made of soft, kink-resistant, NON-DEHP PVC Tubing.  
All products and packaging are not made with natural rubber latex.

Catalog/Part Number	Product Description	
<b>RYM – 7107HPB</b>	<b>InVision-Plus CS® w/Chlorhexidine+Silver Engineering 7” High Pressure Microbore Bifurcated Cath. Ext. Set*</b> <b>(2) Pre-attached InVision-Plus CS®</b> *High Pressure Catheter Extension Set rated 10mL/sec w/max pressure of 325 psi (2) Slide clamps, rotating male-luer lock Note: InVision-Plus CS® is not bonded to set APV = 0.45mL	
<b>RYM – 7307HP</b>	<b>InVision-Plus CS® w/Chlorhexidine+Silver Engineering 7” High Pressure Macrobre Catheter Extension Set*</b> <b>(1) Pre-attached InVision-Plus CS®</b> *High Pressure Catheter Extension Set rated 10mL/sec w/max pressure of 325 psi Slide clamp, universal rotating male-luer lock Note: InVision-Plus CS® is not bonded to set APV = 0.42mL	
<b>RYM – 8104BCL</b>	<b>4” Microbore Bifurcated Catheter Extension Set (2) Pre-attached InVision-Plus® Junior® Clear</b> (2) Slide clamps, rotating male-luer lock Note: InVision-Plus® Junior® is not bonded to set <b>Intended for Use on Neonatal &amp; Pediatric Patients</b> APV = 0.40mL	
<b>RYM – 8107CL</b>	<b>7” Microbore Catheter Extension Set (1) Pre-attached InVision-Plus® Junior® Clear</b> Slide clamp, fixed male-luer lock Note: InVision-Plus® Junior® is not bonded to set <b>Intended for Use on Neonatal &amp; Pediatric Patients</b> APV = 0.27mL	

## Catheter Extension Sets without I.V. Connector

NOTE: All catheter extension sets are made of soft, kink-resistant, NON-DEHP PVC Tubing.  
All products and packaging are not made with natural rubber latex.

Catalog/Part Number	Product Description	
<b>RYM – 3036</b>	<b>36" Ultra-Microbore/Mini Volume Syringe Pump Set</b> Slide clamp, universal rotating male-luer lock APV = 0.30mL	
<b>RYM – 3060</b>	<b>60" Ultra-Microbore/Mini Volume Syringe Pump Set</b> Slide clamp, fixed male-luer lock APV = 0.30mL	
<b>RYM – 3060F</b>	<b>60" Ultra-Microbore/Mini Volume Syringe Pump Set With (1) 0.22 Micron Filter</b> Slide clamp, universal rotating male-luer lock APV = 0.53mL	
<b>RYM – 3104TRF</b>	<b>4" Trifurcated Microbore Catheter Extension Set With (1) 0.22 Micron Filter on 8" tubing segment</b> (3) Slide clamps, universal rotating male-luer lock APV = 0.78mL	
<b>RYM – 3107</b>	<b>7" Microbore Catheter Extension Set</b> Slide clamp, fixed male-luer lock APV = 0.27mL	
<b>RYM – 3360</b>	<b>60" Macrobore Catheter Extension Set</b> Slide clamp, fixed male-luer lock APV = 2.05mL	

Packaging: 100 Units per Box - 4 Boxes per Case - TOTAL: 400 UNITS PER CASE

## Stopcocks

---

NOTE: All catheter extension sets are made of soft, kink-resistant, NON-DEHP PVC Tubing.  
All products and packaging are not made with natural rubber latex.


---

Catalog/Part Number	Product Description
---------------------	---------------------

---

**RYM – 5942-28Y**

**Double 4-Way Stopcock Bonded  
InVision-Plus** with 28" Standard Bore Tubing  
And Y-Injection Site with Check Valve  
APV = 5.75mL



Packaging: 50 Units per Box - 4 Boxes per Case - TOTAL: 200 UNITS PER CASE

---



## ORDERS AND CUSTOMER SERVICE

For Ordering and Customer Service:

Phone: (512) 301-1949  
Fax: (512) 301-7338  
E-mail: customerservice@rymedtech.com

Mail: RyMed Technologies, LLC.  
6000 W. Wm. Cannon Drive  
Building B, Suite 300  
Austin, TX 78749

The minimum order for any product is full case quantities. All orders should be directed to RyMed Technologies LLC Customer Service Department. Customer service hours are 8:00am to 5:00pm, Central Time, Monday through Friday.

### SHIPPING:

Products are shipped domestically via Truck Carrier, Federal Express, or United Parcel Service from Austin, TX or Phoenix, AZ (freight cost borne by RyMed). Express service and/or Drop Shipments are available at the request of the customer (freight cost borne by customer).

### PAYMENT TERMS:

All orders will be billed at current prices at the time of the shipment. Terms are net 30 days for domestic accounts. A service charge of 1 1/2 % per month will be charged to overdue accounts. All prices are in U.S. dollars, and are subject to tax, if any. All prices are subject to change without prior notice. Outside the United States, terms are prepayment by wired funds or letter of credit on sight.

### RETURN GOODS POLICY:

Requests for returns must be addressed to the Customer Service Department. If the reason for return is approved, RyMed Technologies will assign a Return Goods Authorization Number (RGA#). RyMed Technologies will not accept any return without a RGA #. The RGA# must appear on the package returned to RyMed Technologies.

### RETURN GOODS ALLOWANCE POLICY:

- A. Products shipped in error: Full credit, Freight Collect.
- B. Products less than six months old: 75% credit, Dealer pays freight.
- C. Product more than six months old: No credit.
- D. Obsolete or discontinued products (not listed in current price list) - No credit.
- E. Unsalable merchandise returned will be destroyed - No credit.

## OUR MISSION

RyMed is committed to improving patient outcomes aimed at *intraluminal catheter-related bloodstream infection* prevention and *intraluminal thrombotic catheter occlusion* reduction by providing innovative product solutions possessing bona fide fail-safe engineering design technologies.

RyMed's revolutionary products empower healthcare providers to successfully apply proven infection control, catheter care and maintenance practices which support the national goal to deliver healthcare that is safe, cost-effective, and easily adopted.



InVision-Plus®, InVision-Plus CS®, RED®, and Junior® are registered trademarks of RyMed Technologies, LLC.

U.S. Patent Numbers: #6,113,068; 6,299,131; 6,994,315; 7,530,546 B2; 8,096,525, B2 and other U.S. and foreign patents pending.

© 2017 RyMed Technologies, LLC.